Infusomat[®] Space Infusion System Pump Sets

Safety, flexibility and convenience for today's busy care environments



Providing an array of options to encompass your infusion therapy requirements





Infusomat® Space Pump Sets with ULTRASITE® Needle-Free IV System

Description	Product Code	REF Number	Units Per cs
Infusomat Space Pump IV Administration Sets			
		þ	
(15 drops/mL) Universal spike, free flow protection clamp, roller clamp, ULTRASITE injection site 6 inches from distal end, SPIN-LOCK connector. DEHP-free, Latex-free. Priming volume: 23 mL, Length: 120 in. (30)	490031 5 cm)	490031	24
(15 drops/mL) Universal spike, slide clamp, pressure limited check valve, three ULTRASITE injection sites 104 inches, 70 inches and 6 inches above distal end, free flow protection clamp, backcheck valve, roller cla slide clamp, SPIN-LOCK connector. DEHP-free, Latex-free. Priming volume: 22 mL, Length: 120 in. (305 cm)	490032 amp,	490032	24
(15 drops/mL) Universal spike, slide clamp, pressure limited check valve, two ULTRASITE injection sites 88 inches and 6 inches above distal end, free flow protection clamp, backcheck valve at distal Y-site, roller clamp, slide clamp, SPIN-LOCK connector. DEHP-free, Latex-free. Priming volume: 22 mL, Length: 120	490036 in. (305 cm)	490036	24
(60 drops/mL) Universal spike, slide clamp, pressure limited check valve, two ULTRASITE injection sites 82 inches and 6 inches above distal end, free flow protection clamp, roller clamp, flow restrictor, backcheck valve at distal Y-site, SPIN-LOCK connector. DEHP-free, Latex-free. Priming volume: 21 mL, Lengtl	490042 n: 115 in. (292 cm	490042)	24
Infusomat Space Pump Filtered Primary Additive IV Sets			
		e (B)	
(15 drops/mL) Universal spike, slide clamp, pressure limited check valve, two ULTRASITE injection sites 85 inches and 6 inches above distal end, free flow protection clamp, roller clamp, SUPOR® 1.2 micron air-eliminating filter, slide clamp, SPIN-LOCK connector. DEHP-free, Latex-free. Priming volume: 22 mL, Le	490033 ength: 120 in. (30	490033 95 cm)	24
			þ
(15 drops/mL) Universal spike, slide clamp, pressure limited check valve, two ULTRASITE injection sites 82 inches and 6 inches above distal end, free flow protection clamp, roller clamp, SUPOR 0.2 micron air-eliminating filter, slide clamp, SPIN-LOCK connector. DEHP-free, Latex-free. Priming volume: 19 mL, Le	490041 ength: 120 in. (30	490041 95 cm)	24
Infusomat Space Pump Burette Set without Automatic Shutoff			
(60 drops/mL) Non-vented spike, roller clamp, 150 mL burette with 15 micron air guard filter and ULTRASITE injection site, slide clamp, pressure limited check valve, two ULTRASITE injection sites 71 inches and 6 inches above distal end, free flow protection clamp, roller clamp, flow restrictor, backchec at distal Y-site, slide clamp, SPIN-LOCK connector. DEHP-free, Latex-free. Priming volume: 19 mL, Length:		490034	20
Infusomat Space Pump IV Administration Set with ULTRABLOCK [™] UV Resis	tant Tubing		
(15 drops/mL) Universal spike, free flow protection clamp, roller clamp, SUPOR® 0.2 micron air-eliminating	490040	490040	24

(15 drops/mL) Universal spike, free flow protection clamp, roller clamp, SUPOR® 0.2 micron air-eliminating 490040 490040 24 filter, slide clamp, ULTRASITE® injection site 6 inches above distal end,, SPIN-LOCK connector. DEHP-free, Latex-free. Priming volume: 23 mL, Length: 119 in. (302 cm)

Infusomat® Space Pump Sets with ULTRASITE® Needle-Free IV System





Two non-vented spikes, three roller clamps,170 micron blood filter, free flow protection clamp, SAFELINE362043362043injection site six inches above distal end, SPIN-LOCK connector. Latex-free.Priming volume: 39 mL, Length: 120 in. (305 cm)362043

Schedule a demo today.

Call your B. Braun Representative at **1-800-227-2862**

FSC Mixed Sources Product group from weB-managed product gr



www.bbraunusa.com

Supor is a registered trademark of Pall Corporation. Rx Only. ©2010 B. Braun Medical Inc., Bethlehem, PA. All Rights Reserved.